



let's  
keep  
connected



REV 01 02/2019

 polydentia  
swiss manufacture 

Polydentia SA  
via Cantonale 47  
6805 Mezzovico-Vira  
Switzerland

phone: +41 91 946 29 48  
fax: +41 91 946 32 03

info@polydentia.ch  
www.polydentia.ch



# my junior kit

the first  
sectional  
matrix system  
for paediatric  
dentistry

sectional  
matrix  
systems



 polydentia  
swiss manufacture 



new!  
the first  
sectional matrix system  
for paediatric dentistry

# myjunior kit

myJunior kit is the first sectional matrix system for paediatric dentistry, especially developed to suit smaller primary teeth and young permanent teeth. It enables dentists to achieve predictable restorations and ideal contact surfaces. The kit offers child-friendly solutions that help saving chair-time and reduce the distress of children, allowing them to quickly return to their favourite activities.




**myTines Junior**  
Autoclavable and replaceable plastic extremities that fit over the wedge and provide excellent matrix adaptation



**Silicone Tubes**  
Soft extremities that provide extra grip and adaptation



**myWedge**  
Hollow V-shaped, high-tech plastic wedge



**Junior Quickmat Steel Matrices**

use only with rubber dam

kit

**myJunior Kit - REF 6330**

- 1 pc. myClip Junior
- 1 pc. myRing Junior
- 5 pairs Cylindrical Silicone Tubes
- 5 pairs myTines Junior
- 40 pcs. LumiContrast Sectional Matrices, Premolar 0.04mm / 5mm
- 40 pcs. Junior LumiContrast Sectional Matrices, 0.04mm / 3.75mm
- 50 pcs. Quickmat Sectional Matrices, Premolar 0.04mm / 5mm
- 50 pcs. Junior Quickmat Sectional Matrices, 0.04mm / 3.75mm
- 50 pcs. Pinky Wood Wedges, 11mm (XS)
- 50 pcs. Wood Wedges, 12mm (S)
- 50 pcs. myWedges XS
- 50 pcs. myWedges S


refill

- |  |   |   |   |
|--|---|---|---|
| REF 6335 myClip Junior, 2 pcs.   | REF 6336 myRing Junior, 2 pcs.  | REF 6333 myTines Junior, 15 pairs                                     | REF 5711 Silicone Tubes, 15 pairs                                   |
| REF 6708 Junior LumiContrast Sectional Matrices, 0.04mm / 3.75mm 70 pcs. | REF 6705 Junior LumiContrast Sectional Matrices, 0.04mm / 5mm 70 pcs. | REF 5708 Junior Quickmat Sectional Matrices, 0.04mm / 3.75mm 100 pcs. | REF 5705 Junior Quickmat Sectional Matrices, 0.04mm / 5mm, 100 pcs. |
| REF 5011 Pinky Wood Wedges, 11 mm (XS), 200 pcs.                         | REF 5012 Wood Wedges, 12 mm (S), 100 pcs.                             | REF 7101 myWedge, Extrasmall (XS), 100 pcs.                           | REF 7102 myWedge, Small (S), 100 pcs.                               |

myClip Junior and myRing Junior are paediatric sectional matrix rings that feature minimally invasive, safe to use, bright-coloured solutions that allow easier restorative treatment on young teeth of low crown height. The junior sectional matrix rings combine stability with the right separation force enabling dentists to achieve ideal contact surfaces. In addition, the extra grip provided by the flexible extremities myTines Junior and rounded silicone tubes ensure a better adaptation of the sectional matrices, that helps preventing annoying composite flashes and overhangs and saving valuable chair-time.

## myClipjunior

All in one paediatric sectional matrix ring with integrated forceps



ADVANTAGES

- All in one sectional matrix ring. No need for intimidating forceps
- Long-lasting ring, made out of Nickel Titanium
- Optimal separation force to achieve flawless contact surfaces
- Quick, easy, one-hand ring positioning
- Autoclavable and replaceable plastic myTines Junior extremities



"I regard myClip Junior as an optimal solution for paediatric treatments as it allowed me to achieve an ideal restoration with less efforts and with considerable chair-time savings."


Dr. Marina Papachroni  
Greece



for more information


## myRingjunior

Paediatric sectional matrix ring that fits a wide variety of clinical situations



ADVANTAGES

- Ideal size to fit on deciduous teeth in children's mouth
- Long-lasting ring made out of Nickel Titanium
- Optimal separation force to achieve flawless contact surfaces
- Excellent matrix adaptation and retention
- Autoclavable and replaceable plastic myTines Junior extremities



"Children restorative treatments won't be tricky business anymore thanks to myRing Junior and the paediatric solutions from Polydentia. Paediatric restorations are easy and efficient thanks to the small dimension of the sectional matrix ring and the junior matrices. The system ensures a perfect fit on milk teeth and prevent time-consuming composite flashes. Restorations are also more enjoyable as my little patients don't have to withstand prolonged restoration times."

Dr. Martina Sabová  
Slovakia



for more information